

# LINC Night Bags



## LINC-Flo Reusable

LINC-Flo reusable night drainage bags allow problem free drainage to help patients have a good nights' rest.



Long flexible tubing and needle free sample port

Integrated bed hanger

**TZ02**

Approx vol. ml.

2000  
1800  
1600  
1400  
1200  
1000  
800  
600  
400  
300  
200  
100

100  
75  
50  
25

GLOVES

T-tap and lever tap available

**SLC13**



Each bag includes:



Nitrile gloves

## Single use Night Bags



LINC Night bags come with either a single use or reusable taps.

## Instructions for Use:

### Connection

- Wash hands. Take out the bag from its packaging and check that the taps are in the closed position.
- Remove connector cap and insert the ridged connector fully for a secure connection into the catheter funnel for reusable taps or leg bag funnel for single use taps.
- Ensure that the night bag is positioned below the leg bag or bladder and supported by an appropriate hanger or stand
- Ensure that any of the tubing is NOT kinked and the leg bag tap is open (if connected to one).
- All LINC Night Bags are for single use only and must not be reused once disconnected.

### Emptying

- Empty the night bag by opening the tap over a toilet or designated receptacle.
- All LINC bags are ready to use and come with the taps closed
- **Lever tap:** open the tap by pushing the lever downwards towards the outlet tap and close by pulling the lever tap upwards towards the bag
  - o The lever points towards the desired flow of urine
  - o Tap starts to open/close after the 90° angle of travel
  - o The whole tap section may be twisted to avoid accidental opening
- **T-tap:** open and close the tap by pushing the barrel sideways
- Remember to close the tap after emptying.
- Ensure that the outlet tap is as dry and without any residual urine as possible
- **Twist tap:** turn the bag upside down with the tap pointing up; using both hands, hold the tap and twist-off the outermost tip; invert the bag to allow the urine to flow out into the toilet/receptacle

### Usage and Disposal

- LINC night bags with re-usable taps should be changed between 5 to 7 days of use.
- LINC night bags with single use taps should be disposed of after emptying.
- Store at room temperature
- Dispose with normal household waste or as directed by your clinician
- Do not flush down the toilet

Tap	type	Use	Sterile	Capacity
Twist Tap	single use	non sterile	2L	
	single use	sterile	2L	
	single use	sterile	3L	
T-Tap	single use	sterile	2L	
	re-usable	sterile	2l	
	re-usable	sterile	3L	
Lever Tap	re-usable	sterile	2L	

**Distributed** in the UK by LINC Medical. **Manufactured by** Changshu Wilson Medical Devices Co., Ltd., No 1 Taikang road, North of Zhitang Town, Changshu City, 215531 Jiangsu Province, China. **EC Representative:** Welkang Ltd, Suite B, 29 Harley Street, London W1G 9QR, UK.