# **Continence Products** Community order codes




Leading the way in **continence care** 

### Who are LINC?

LINC Medical is a UK based, family owned continence product supplier. Because of our size and understanding of UK healthcare needs, we rapidly innovate to bring life changing products to patients and healthcare professionals.

## Contents

### Co CATHETERS



### DRAINAGE BAGS

12 LINC-Flo Leg Bag

ľ

- 13 LINC-Flo Nephrostomy Leg Bag
- 14 Panda Bag Paediatric Leg Bag
- 15 LINC Leg Bag Sleeve
- 16 LINC-Flo Reusable Night Bags
- **17** LINC-Flo Single-Use Night Bags

### ACCESSORIES

18 Fixstrip

*C*JP

- 18 Care-Flo Catheter Valve
- **19** Bladder Infusion Kit (BIK)
- 20 Common Questions
- 21 Catheter Troubleshooting

# Why use a LINC catheter?



#### Specialist Catheter Tip Designs

Patients experience different catheterrelated problems such as blockages, spasms and recurrent infection. Our unique catheter designs allow you to apply individualised patient care, choosing the right catheter tip to suit the patient.



#### Integral Balloon Technology

The LINC catheter balloon is integral to the wall of the catheter. This reduces the pain and trauma associated with balloon cuffing upon removal which can occur with standard catheter balloons.<sup>2,4</sup>





12ch

0

#### 100% Medical Grade, Biocompatible Silicone (Latex free)

Schematic Integral Balloon

Silicone catheters provide double the drainage compared to the equivalent sized latex catheters.<sup>4,6</sup> This helps to prevent blockages and ensure excellent, continuous drainage.<sup>4,7</sup>

Ō



#### Scientifically Validated Inflation Solution

LINC catheters are supplied with a pre-filled syringe of 5% glycerine solution that is proven to reduce premature balloon deflation. This solution equalises the osmotic pressure between the balloon and urine and helps the catheter to stay in situ for the maximum recommended time.<sup>4,8</sup>



Optitip

### Small changes make big differences

The Optitip catheter has been designed with a balloon that is closer to the tip. This prevents damage to the bladder mucosa that can result from standard catheter tips which may suck against the bladder wall. An eyelet has been added below the balloon to drain residual urine.

Hole below balloon encourages drainage of residual urine

•

Healthy bladder mucosa

> Residual urine drainage

Includes syringe of Gel Lubricant (with lidocaine and chlorhexidine) Balloon moved closer to tip to reduce bladder wall damage

> Excellent drainage

centro

) 100% Silicone 🕢 Integral Balloon 🕢 Two Wa

| 6

"I no longer have bladder pains from pointed catheter tips and therefore have no catheter induced spasms of bladder, legs or bowels."



**Patient** Frimley Park, Surrey

### ADULT CODES

СН	Product Code	Description	Balloon	PIP code
<b>O</b> 12CH	1007UB12OT	40cm Standard Length	5ml	4045910
<b>1</b> 4CH	1007UB14OT	40cm Standard Length	5ml	4045928
<b>1</b> 6CH	1007UB16OT	40cm Standard Length	5-10ml	4045936
<b>1</b> 8CH	1007UB18OT	40cm Standard Length	5-10ml	4045944
O 12CH	1007UBF12OT	28cm Female Length	5ml	4045951
<b>1</b> 4CH	1007UBF14OT	28cm Female Length	5ml	4045969
<b>1</b> 6CH	1007UBF16OT	28cm Female Length	5-10ml	4045977
18CH	1007UBF18OT	28cm Female Length	5-10ml	4045985

#### > PAEDIATRIC CODES

СН	Product Code	Description	Balloon	PIP code
<b>8</b> CH	1007UBS08OT	Paediatric 33cm	3ml	4045894
• 10CH	1007UBS10OT	Paediatric 33cm	3ml	4045902

# Open Tip

The open tipped catheter provides excellent drainage with both an open-ended tip and two drainage holes above the balloon.

> Includes **free** syringe of Lubricant (with lidocaine and chlorhexidine)



СН	Product Code	Description	Balloon	PIP code
<b>O</b> 12CH	08451205G	40cm Standard Length	5ml	3834066
<b>1</b> 4CH	08451405G	40cm Standard Length	5ml	3834074
<b>1</b> 6CH	08451610G	40cm Standard Length	5-10ml	3834082
<b>1</b> 8CH	08451810G	40cm Standard Length	5-10ml	3834090
<mark>0</mark> 20CH	08452010G	40cm Standard Length	5-10ml	3834108
<b>2</b> 2CH	08452210G	40cm Standard Length	5-10ml	3834116
<b>2</b> 4CH	08452410G	40cm Standard Length	5-10ml	3834124
O 12CH	08471205G	28cm Female Length	5ml	3834132
14CH	08471405G	28cm Female Length	5ml	3834140
<b>1</b> 6CH	08471610G	28cm Female Length	5-10ml	3834157
18CH	08471810G	28cm Female Length	5-10ml	3834165
<b>20CH</b>	08472010G	28cm Female Length	5-10ml	3834173



СН		Product Code	Description	Balloon	PIP code
	8CH	08470803G	Paediatric 28cm	3ml	3833910
•	10CH	08471003G	Paediatric 28cm	3ml	3833928

| 8

11º

 $\bigcirc$ 

# Tiemann Tip

The tiemann tip catheter has a unique shape that allows easier insertion in patients who have an enlarged prostate or urethral stricture.



 $\langle \rangle$ 

Includes **free** syringe of Lubricant (with lidocaine and chlorhexidine)

Two Way

100% Silicone 🕢 Integral Balloon 🔗

## ADULT CODES

СН	Product Cod	e Description	Balloon	PIP code
<b>O</b> 12CH	1011UB12	40cm Standard Length	5ml	3945193
<b>1</b> 4CH	1011UB14	40cm Standard Length	5ml	3886777
<b>1</b> 6CH	1011UB16	40cm Standard Length	5-10ml	3886751
<b>1</b> 8CH	1011UB18	40cm Standard Length	5-10ml	3886769



 $\bigcirc$ 

# **Round Tip**



Includes **free** syringe of Lubricant (with lidocaine and chlorhexidine)

100% Silicone 🕢 Integral Balloon\* 📿 Two W

Feature not applicable to 6CH paediatric Round Tip catheters

ADULT CODES

СН	Product Code	Description	Balloon	PIP code
<b>O</b> 12CH	08501205G	40cm Standard Length	5ml	3834280
<b>1</b> 4CH	08501405G	40cm Standard Length	5ml	3834298
• 14CH	08501410-OTG	40cm Standard Length	10ml	4137311
<b>—</b> 16CH	08501610G	40cm Standard Length	5-10ml	3834306
<b>1</b> 8CH	08501810G	40cm Standard Length	5-10ml	3834314
🛑 12CH	085012051G	33cm Female Length	5ml	3834322
• 14CH	085014051G	33cm Female Length	5ml	3834330
14CH	08501410F-OTG	33cm Female Length	10ml	4137329
○ 16CH	085016101G	33cm Female Length	5-10ml	3834348
<b>O</b> 14CH	08501430	40cm Standard Length	20ml	2828325
<b>1</b> 6CH	08501630	40cm Standard Length	30ml	2828333

### PAEDIATRIC CODES

Cŀ	H P	roduct Code	Description	Balloon	PIP code	Stylet
	6CH	08500601G	Paediatric 33cm	1.5ml	3834256	Non-Stylet
	6CH	08500601SG	Paediatric 33cm	1.5ml	4012472	Stylet
	8CH	08500803G	Paediatric 33cm	3ml	3834264	Non-Stylet
	8CH	08500803SG	Paediatric 33cm	3ml	4012480	Stylet
	10CH	08501003G	Paediatric 33cm	3ml	3834272	Non-Stylet
	10CH	08501003SG	Paediatric 33cm	3ml	4012498	Stylet

# Round Tip

The LINC Round Tip Xtra utilises our standard round tip technology for easy insertion and adds a drainage hole below the balloon to aid residual urine drainage.



 $\langle \rangle$ 

Includes **free** syringe of lubricant (with lidocaine and chlorhexidine)

Two Way

100% Silicone 🕢 Integral Balloon 🗸

### ADULT CODES

СН	Product Code	Description	Balloon	PIP code
<b>O</b> 12CH	1007UB12	40cm Standard Length	5ml	4046637
<b>1</b> 4CH	1007UB14	40cm Standard Length	5ml	4046645
<b>1</b> 6CH	1007UB16	40cm Standard Length	5-10ml	4046652
<b>1</b> 8CH	1007UB18	40cm Standard Length	5-10ml	4046660
O 12CH	1007UBF12	28cm Female Length	5ml	4046678
14CH	1007UBF14	28cm Female Length	5ml	4046686
<b>1</b> 6CH	1007UBF16	28cm Female Length	5-10ml	4046694
18CH	1007UBF18	28cm Female Length	5-10ml	4046702

## PAEDIATRIC CODES

СН		Product Code	Description	Balloon	PIP code
	8CH	1007UBS08	Paediatric 33cm	3ml	4046710
	10CH	1007UBS10	Paediatric 33cm	3ml	4046728

-**Line**-

## LINC-Flo Leg Bag

LINC-Flo sterile leg bags offer proven patient comfort, privacy & safety.

> Soft fabric backing

Easy one-handed operation of tap

Needle free port protects against infection

Adjustable bag includes:

ach bag ncludes:

80

300

20.

0.

50

Latex free gloves Scissors

Straps

12



Product Code	Capacity	Tube	Тар	Box Qty	PIP code
LM350SD-T	350ml	Direct	Т	10	3384997
LM350MD-T	350ml	10cm	Т	10	3385002
LM500SD-T	500ml	Direct	Т	10	3383726
LM500MD-T	500ml	10cm	Т	10	3385036
LM500LD-T	500ml	30cm	Т	10	3385044
LM500AD-T	500ml	ADJ	Т	10	3385051
LM500SD-L	500ml	Direct	Lever	10	3712114
LM500MD-L	500ml	10cm	Lever	10	3712213
LM500LD-L	500ml	30cm	Lever	10	3783891
LM500AD-L	500ml	ADJ	Lever	10	3811643
LM750SD-T	750ml	Direct	Т	10	3383692
LM750MD-T	750ml	10cm	Т	10	3383700
LM750LD-T	750ml	30cm	Т	10	3383718
LM750AD-T	750ml	ADJ	Т	10	3385085
LM750AD-L	750ml	ADJ	Lever	10	3834363



### LINC-Flo Nephrostomy Leg Bag

A comfortable, fabric backed leg bag with a nephrostomy connector and pair of latex free gloves

Product Code	LM750NDB
Capacity	750ml
Тар	Т
Box Qty	10
PIP code	4044103

--

## Panda Bag Paediatric Leg Bag

Panda bags cater to the needs of infants and babies who require urine drainage.

Pleasing, child-friendly panda-themed design

- Available with direct inlet, 10 cm tubes, 30 cm tube and in 150 ml and 250 ml capacities to suit different ages and needs.
- Soft fabric backing ensures the panda bag is comfortable
- Needle-free sample port eliminates risk of needlestick injury
- Each bag is supplied with an elasticated strap and pair of latex free gloves



Product Code	Capacity	Tube	Тар	Box Qty	PIP code
LM100SD-L	100ml	Direct	Lever	10	4050969
LM100MD-L	100ml	10cm	Lever	10	4050993
LM100LD-L	100ml	30cm	Lever	10	4051009
LM200SD-L	200ml	Direct	Lever	10	4051017
LM200MD-L	200ml	10cm	Lever	10	4050977
LM200LD-L	200ml	30ml	Lever	10	4050985



 $\bigcirc$ 

## LINC Leg Bag Sleeve

A comfortable, discreet and secure alternative to straps.

### ADULT CODES

Product Code	Description	Qty/pack	PIP code
ME350 750S S	SMALL (leg circum 24-39cr	n) 4	4099552
ME350 750S M	MEDIUM (leg circum 36-55	cm) 4	4099537
ME350 750S L	LARGE (leg circum 40-70cr	n) 4	4099529
ME350 750S XL	X-LARGE (leg circum 65-95	cm) 4	4099511
ME350 750S XXL	XX-LARGE (lea circum 75-1	05cm) 4	4099750



## LINC-Flo Reusable

LINC-Flo reusable night drainage bags allow problem free drainage to help patients have a good nights' rest.

Long flexible tubing and needle free sample port Integrated bed hander T-tap and lever tap available 2000 1800 1600 1400 1200 1000 Approx vol. ml. 800 600 Room SLC13 Detun/date:

**Product Code** 

TZ02

TZ02-L

SLC13

SLC13-L

Capacity Tap

Т

Т

Lever

Lever

2L

2L

2L

2L

100

**Box Oty** 

10

10

10

10

**PIP code** 

3054822

3720455

3406071

3712247

TZ02
Vent
Non return va

Each bag

includes:


Latex free

gloves

SLC13

--





### LINC-Flo Single Use 2L and 3L

The LINC-Flo single use night drainage bag provides a clean, safe and easy solution for night urine drainage.

Product Code	LM2LNS	LM2LS	LM3LS
Description	Non-sterile	Sterile	Sterile
Capacity	2L	2L	3L
Box Qty	10	10	10
PIP code	3054780	3054806	3945201

# LINC-Flo Single Use 2L and 3L with T-tap

The T-tap allows for easier emptying and gives patients' choice.

Product Code	LM2LS-T	LM3LS-T
Description	Sterile	Sterile
Capacity	2L	3L
Box Qty	30	30
PIP code	4131744	4131736



### **FixStrip**

A customisable catheter fixation strap

100cm length allows flexibility for patient to choose whether it is applied to the thigh or abdomen

) Free pair of scissors to cut FixStrip to desired length

Product	Box	PIP
Code	Qty	code
FIXSTRIP	4	4049938

### Care-Flo Catheter Valve

A 'lever tap' valve which securely fits all types of indwelling catheters.

Handy lever mechanism makes it easy for patients with poor manual dexterity to initiate voiding.

Firm mechanism helps to prevent accidental valve opening

Ridged connector ensures a tight fit with all brands of Foley Catheter

Product	Box	PIP
Code	Qty	code
CF1	5	3441524

### (BIK) Bladder Infusion Kit

Allows clinicians to easily administer catheter maintenance solution without the need to disconnect the catheter and drainage bag. Single use.

Aseptic technique can be maintained and infection risk can be minimised

Clamp for stopping flow to the leg bag

Funnel which provides access to needle-free sample port

Includes sterile wipes to clean the sample port before use



Product	Box	PIP
Code	Qty	code
MCI/701	10	

V

# Your catheter and your lifestyle: common questions

#### Can I go swimming with my catheter?

You can go swimming with either a urethral or suprapubic catheter. Remember in the case of a suprapubic catheter a small incision will have been made, so swimming should be kept until at least a few days after the operation. Make sure the skin around the catheter site has fully healed and shows no signs of infection before getting in the water.

#### Can I have sex with a catheter in place?

In the case of suprapubic catheters, this is in fact one of the advantages compared to a urethral one in that your genitalia are kept free for sex. Sex is still possible with a urethral catheter but may be more of a challenge.

Women can fold the catheter up to one side of the abdomen and tape it down during sex. Men can fold the catheter along the side of the penis and use a condom to keep in place. If you know how to remove and replace a catheter then this may also be an option after discussion with your doctor or nurse.

Before having sex partners should wash to reduce the chance of infection. To reduce the risk of damage to either the urethra or site on the abdomen rough sex should be avoided.

It is very important to not use a lubricant such as Vaseline

or any other petroleum jelly based lubricant. These lubricants can damage the catheter tubing. Use lots of water-soluble lubricant instead like KY jelly.



### What should I do if...

#### CATHETER TROUBLESHOOTING

# 1. My patient has bladder spasms

Spasms can be caused by the pointed catheter tip irritating the bladder wall, an over inflated balloon or over-sized catheter diameter irritating the bladder trigone and over concentrated urine which may irritate the bladder wall.<sup>9</sup>

- Ensure the patient is well hydrated
- Try a smaller sized catheter and/or a catheter with a smaller balloon
- Trying using an Optitip Catheter to reduce bladder wall irritation
- Check that the bag is not pulling on the catheter. Consider using a fixation device such as FixStrip or LockBand.

#### 2. My patient experiences trauma and pain upon insertion/removal

Trauma and pain are often caused by cuffing of the catheter balloon once deflated.<sup>2</sup> These issues can make catheter removal problematic in patients who have a suprapubic or urethral catheter.

- Try using a catheter with an integral balloon, which maintains the same diameter of the catheter along the whole length and prevents cuff formation.<sup>3,4</sup>
- Ensure the smallest Charriere catheter possible is used for the patient (not applicable to suprapubic insertion).

# **3.** My patient consistently has poorly draining catheters

Poor drainage may be caused by mechanical problems such as kinked tubing and bag positioning, or luminal blockage due to encrustation of the catheter.<sup>10</sup>

- Ensure there are no kinks in the catheter or bag tubing
- Ensure the drainage bag is below the bladder
- Try using a 100% silicone catheter which has a wider drainage lumen and better Flows'
- Try using a catheter that is designed For high drainage, such as the LINC Open Tip.

### What should I do if...

#### **CATHETER TROUBLESHOOTING**

# **4.** My patient is bypassing

Bypassing may result from mechanical blockage (problem 3) and irritation of the bladder trigone (problem 1).

Another reason may be the kinking of the catheter around the eyelets due to abdominal pressure.<sup>9</sup>

- Try the Optitip catheter which has minimal, open tipped design to prevent kinking at the eyelets.
- Follow the advice from problems 1 and 3.

# **5.** My patient is suffering from recurring UTI's

Having a long-term urinary catheter increases the risk of developing urinary tract infections (UTIs). There may be several causes of UTIs such as poor hygiene, a weakened bladder wall due to catheter tip irritation or an infected sump of residual urine.

- Ensure the patient is following correct hygiene practice, remains hydrated, consumes a high fibre diet and regularly checks there are no kinks in the catheter.
- Consider switching the patient to a catheter that helps to drain residual urine and protects the bladder wall from damage such as the Optitip.

# **6.** My Patient's catheter regularly falls out

Catheters may fall out when balloons gradually deflate as fluid leaves the balloon due to differences in osmotic pressure between the inflation solution and the urine.<sup>4</sup>

 Try using an Osmotically balanced inflation solution such as the 5% glycerine solution supplied with LINC Medical Catheters which helps to maintain balloon inflation for the duration of catheterisation.<sup>11</sup>

#### References

Robinson, J. Suprapubic catheter removal: the cuffing effect of deflated catheter balloons. Br. J. Community Nurs. 8, 205–208 (2003).
 Leaver, R. Understanding long-term catheterisation for effective bladder drainage. J. Community Nurs. 31, 43–50 (2017).
 Robinson, J. 100% Silicone Catheters with itegral balloon for suprapubic catheterisation. (2005).
 Lawrence, E. L. & Turner, I. G. Kink, flow and retention properties of urinary catheters part 1: Conventional foley catheters. J. Mater. Sci. Mater. Med. 17, 147–152 (2006).
 Robinson, J. Selecting a urinary catheter and Drainage System. Br. J. Nurs. 15, 1045–1049 (2006).
 Barford, J. M. T., Hu, Y., Anson, K. & Coates, A. R. M. A Biphasic Response From Bladder Epithelial Cells induced by Catheter Material and Bacteria: An In Vitro Study of the Pathophysiology of Catheter Related Urinary Tract Infection. J. Urol. 180, 1522–1526 (2008).
 Luck M. Catheter usage: a design challenge. Urol. Nurs. 32, 307–312 (2012).
 Evans, A., Painter, D. & Feneley, R. Blocked urinary catheters: nurses' preventive role. Nurs. Times 97, 37 (2001).
 Aslam, N. & Moran, P. A. Catheter use in gynaecological practice. Obstet. Gynaecol. 16, 161–168 (2014).



Family run



Rapid innovation



British Company



Address: LINC Medical Systems Lto 7 Kingsley Street, Leicester, LE2 6DY

**Telephone:** +44 (0) 1572 717515 **Email:** sales@linc-medical.co.ul

Website: www.linc-medical.co.uk

© 2018 LINC Medical CMB V4 270922