

GENESIS™

TUBE SEALERS



RAPID SEAL II™
SE640
PORTABLE UNIT



RAPID SEAL II™
SE530
MULTI-HEAD SEALER



RAPID SEAL II™
SE540
SPACE SAVING BENCHTOP SEALER



RAPID SEAL II™
SE330 & SE340
STANDARD BENCHTOP UNIT

Component processing is too important to leave to chance.

SPEED

- No warm-up required – Instant on
- Rapid Seal II Technology
- Never overheat
- Notched Seal™ for easy separation

COMPATIBILITY & RELIABILITY

- “Intelligent” RF Board – accepts and adjusts to all size blood tubing
- Legendary Genesis brand quality
- Trusted by blood banks worldwide
- 1 year warranty on parts and labor
- CE and UL compliant


VERSATILITY & SAFETY

- Standard bench-top units
- Portable units (AC/Battery)
- Space saving units
- Flexible multi-head solutions

GENESIS™ TUBE SEALERS

WITH RAPID SEAL II TECHNOLOGY



Model	SE330	SE340	SE370	SE530	SE540	SE640	
General Description	STANDARD BENCHTOP SEALER	MULTIPLE STATION BENCHTOP SEALER	BATTERY POWERED BENCHTOP SEALER	SPACE SAVING BENCHTOP SEALER	SPACE SAVING BENCHTOP SEALER	PORTABLE SEALING UNIT	
Sealing Head Configuration	Fixed sealing head	Fixed and Hand-held sealing unit	Fixed and Hand-held sealing unit	Fixed sealing head	Fixed and Hand-held sealing unit	Hand-held sealing unit	
Catalog #	428-SE330	428-SE340	428-SE370	428-SE530	428-SE540	428-SE640	
Special Features	<ul style="list-style-type: none"> Retractable splash guard “Intelligent” RF Board RAPID SEAL II seals in 0.5 to 1.2 seconds Easy Clean head protector 	<ul style="list-style-type: none"> Retractable splash guard 6’ standard cable “Intelligent” RF Board RAPID SEAL II seals in 0.5 to 1.2 seconds Easy Clean head protector Use fixed and hand-held sealing unit independently 	<ul style="list-style-type: none"> Retractable splash guard 6’ standard cable “Intelligent” RF Board RAPID SEAL II seals in 0.5 to 1.2 seconds Battery Operated and rechargeable Easy Clean head protector 2000+ Seals per full battery charge Aluminum carrying case 	<ul style="list-style-type: none"> Flexible multi-head option Splash guard “Intelligent” RF Board RAPID SEAL II seals in 0.5 to 1.2 seconds Custom Sealer Stand Easy Clean head protector Slim design to maximize available space Flexible Head Spacing (75 / 80 / 85 / 90 mm) 	<ul style="list-style-type: none"> Flexible multi-head option Splash guard “Intelligent” RF Board 6’ standard cable RAPID SEAL II seals in 0.5 to 1.2 seconds Custom sealer stand Easy Clean head protector Slim design to maximize available space Flexible Head Spacing (75 / 80 / 85 / 90 mm) 	<ul style="list-style-type: none"> 6’ standard cable RAPID SEAL II seals in 0.5 to 1.2 seconds Hand Sealing Unit Battery operated and rechargeable Lightweight canvas carrying bag included Medical device classified* 1000+ Seals per full battery charge 	
DIMENSIONS	Main Unit <small>Inches mm</small>	6.9W x 6.1H x 13D 176W x 156H x 331D	6.9W x 6.1H x 13D 176W x 156H x 331D	6.9W x 6.1H x 13D 176W x 156H x 331D	2.9W x 9.7H x 12.8D 74W x 248H x 326D	2.9W x 9.7H x 12.8D 74W x 248H x 326D	6.1W x 1.5H x 14.4D 156W x 37H x 265D
	Hand Held Sealing Unit <small>Inches mm</small>	N/A	1.1W x 1.8H x 8.1D 28W x 44H x 204D	1.1W x 1.8H x 8.1D 28W x 44H x 204D	N/A	1.1W x 1.8H x 8.1D 28W x 44H x 204D	1.1W x 1.8H x 8.1D 28W x 44H x 204D
	Unit <small>lbs Kg</small>	11.0 lbs 5.0 Kg	11.0 lbs 5.0 Kg	11.0 lbs 5.0 Kg	9.7 lbs 4.4 Kg	9.7 lbs 4.4 Kg	4.0 lbs 1.8 Kg
	Sealing Gun <small>lbs Kg</small>	N/A	0.5 lbs 0.2 Kg	0.5 lbs 0.2 Kg	N/A	0.5 lbs 0.2 Kg	0.5 lbs 0.2 Kg
Power Type	A/C	A/C	Battery Operated	A/C	A/C	Battery Operated	
Power Source	100-240 V 50/60 Hz	100-240 V 50/60 Hz	12V / 10AH NiMH Battery (rechargeable)	100-240 V 50/60 Hz	100-240 V 50/60 Hz	12V / 4.5AH NiMH Battery (rechargeable)	
Max Tubing Dia.	6 mm (1/4")	6 mm (1/4")	6 mm (1/4")	6 mm (1/4")	6 mm (1/4")	6 mm (1/4")	
Remarks	—	—	—	 <p>The 530 & 540 sealers can be used stand alone or linked as a multi-head sealing unit. Link 3, 4, 5, 6 or more units together (High Volume Use)</p>		*Approved for use during a standard donation or apheresis procedure.	

*All devices are CE and UL compliant.

CORPORATE HEADQUARTERS:
65 Commerce Way
Hackensack, NJ 07601
TEL: 201-708-1400
TEL: 866-71-BLOOD
FAX: 201-708-1104
info@GenesisBPS.com

INTERNATIONAL:
Genesis BPS International Sàrl
CP 2361, 1260 Nyon 2, Switzerland
international@GenesisBPS.com

www.GenesisBPS.com

GENESIS BPS™

Technology is in our Blood.™