

GENESIS[™] BPS[™]

Specialty Cell Processing Products

Safely seal cord blood units using the Genesis Cord Blood Sealer and Cord Blood Sealing Jig.

GENESIS Cord Blood Sealer

CAT# 428-SE470RS

Product Features

- Specially designed sealing head for use with the Pall[™] cord blood freezing-bag chamber
- Prevents cell necrosis when used with Genesis Cord Blood Sealing Jig
- Benchtop space saving design
- Multiple station use (6' cable)
- Carrying bag fits all accessories
- Max tubing diameter is 6mm (1/4")



GENESIS Sealing Jig

CAT# 428-UCBJIG

Product Features

- Designed to work with the Pall[™] 25mL cord blood freezing-bag chamber
- Accurately guides tube sealing head
- Protects freezing-bag corners
- Designed to prevent arching
- Exclusive to GenesisBPS

To Order Contact GenesisBPS

465 Route 17 South, Ramsey, NJ 07446
TEL: 866-71-BLOOD / 201-708-1400 FAX: 201-708-1104
EMAIL: info@GenesisBPS.com

GENESIS[™] BPS[™]

Where the global bloodstream beginsSM

www.GenesisBPS.com